

Name of the medical device

Sensodyne Rapid Relief Mint Toothpaste

Sensodyne Rapid Relief Whitening Toothpaste

Intended purpose

Sensitivity Relief

Directions for use

Brush twice a day and not more than three times, minimise swallowing and spit out. For fast relief, brush sensitive teeth first or squeeze a pea-sized amount onto a clean fingertip and gently rub onto each sensitive area (max twice per day) for 1 minute before brushing. Replace lid to keep product dry.

Warnings and precautions

Keep out of reach of children. Not for use by children under 12 years of age, unless on advice of a dentist or doctor. Sensitive teeth may indicate an underlying problem that needs prompt care, if symptoms persist or worsen see your dentist. If irritation occurs, discontinue use. Do not use if foil seal is broken. Do not use if carton ends are open. Report any product incident or side effects to the manufacturer, your Competent Authority and a pharmacist or dentist.

Manufacturer

Legal manufacturer: Haleon Medical Devices, Clocherane, Youghal Road, Dungarvan, Co. Waterford, X35 Y983, Ireland.

UK Responsible Person: Haleon UK Trading Limited, The Heights, Weybridge, KT13 0NY, U.K.