

✓BLEEDING ✓REDNESS ✓SWELLING

CLINICALLY PROVEN TO HELP REVERSE SIGNS OF EARLY GUM DISEASE:



Plaque buildup on gums and teeth can lead to gum disease. parodontax Active Gum Repair is a daily fluoride toothpaste formulated to help reverse signs of early gum disease by targeting plaque, the main cause of gum problems. It is 3xs more effective at removing plaque*

With a fresh minty flavor to enjoy every day.

*vs a sodium monofluorophosphate toothpaste with regular brushing

Distributed by: Haleon, Warren, NJ 07059 @2024 Haleon or licensor.

Targets Improves gum health

1 1 1 1 1 1

1 1 1 1 1 1

1 1 1 1 1 1

1 1 1 1 1 1





Drug Facts

Active ingredient
Stannous fluoride 0.454% (0.15% w/v fluoride ion).

Purposes

.. Anticavity, Antigingivitis

Uses • aids in the prevention of dental cavities. • helps control bleeding gums. helps interfere with harmful effects of plaque associated with gingivitis.

Warnings

Stop use and ask a dentist if • gingivitis, bleeding, or redness persists for more than 2 weeks.
• you have painful or swollen gums, pus from the gum line, loose teeth, or increasing spacing between the teeth. These may be signs or symptoms of periodontitis, a serious form of gum disease.

Keep out of reach of children under 6 years of age. If more than used for brushing is accidentally swallowed, get medical help or contact a Poison Control Center right away.

Directions • adults and children 2 years of age and older: • brush teeth thoroughly, preferably after each meal at least twice a day, and not more than 3 times a day, or as directed by a

Drug Facts (continued)

dentist or doctor. Minimize swallowing. Spit out after brushing. • to minimize swallowing for children under 6 years of age, use a pea-sized amount and supervise brushing and rinsing until good habits are established. • children under 2 years of age: Consult a dentist or doctor.

Other information • products containing stannous fluoride may produce surface staining of the teeth. Adequate toothbrushing may prevent these stains which are not harmful or permanent and may be removed by your dentist. • this product is specially formulated to help prevent staining. do not store above 25°C (77°F)

 $\textit{Inactive ingredients} \quad \text{glycerin, PEG-8, hydrated silica, pentasodium triphosphate, flavor,}$ sodium lauryl sulfate, titanium dioxide, polyacrylic acid, cocamidopropyl betaine, sodium saccharin

Questions or comments? 1-866-844-2797

ALWAYS FOLLOW THE LABEL

1 1 1 1 1

1 1 1 1 1 1

1 1 1 1 1 1

1 1 1 1 1 1

